

Built on a heritage of proven performance,
the refined SharpPoint Clear Control 3D™ Plus is here!

The Refined SharpPoint™ Clear Control 3D™ Plus Knives Are Now Available.

CC3D™ Plus slit knives use a world-class blade tip with uniform strength and ductility that is critical for creating a smooth, well-controlled incision using less force than a traditional phaco slit knife.

**SHARPNESS + STRENGTH =
IMPROVED PERFORMANCE**

100% MADE IN THE U.S.A.



Sharpness and Control

- 33 % lower penetration force required than standard slit knives¹
- Blade design allows for less vertical tissue displacement for a more self-sealing wound
- Thinner blade decreases distortion of tissue resulting in improved wound architecture
- Excellent control and minimal drag

3D Incision Design


- Multi-dimensional incision
- Tip design adds vertical dimension, placing less stress on wound margins
- Consistent use of force no matter what size CC3D knife

CC3D Plus Slit Knives

ORDERING INFORMATION

ITEM NO.	DESCRIPTION	KNIVES PER BOX
72-1650	CC3D 1.6mm Knife, Angled	6
72-1850	CC3D 1.8mm Knife, Angled	6
72-2050	CC3D 2.0mm Knife, Angled	6
72-2250	CC3D 2.2mm Knife, Angled	6
72-2450	CC3D 2.4mm Knife, Angled	6
72-2550	CC3D 2.5mm Knife, Angled	6
72-2650	CC3D 2.65mm Knife, Angled	6
72-2750	CC3D 2.75mm Knife, Angled	6
72-2850	CC3D 2.85mm Knife, Angled	6
72-3050	CC3D 3.0mm Knife, Angled	6
72-3250	CC3D 3.2mm Knife, Angled	6

PENETRATION STUDY



¹ Comparative study on file of SharpPoint knives. Average values.

CC3D Plus Blade

- *Stealth-like tip and thinner blade prevents tenting during incision resulting in superior wound architecture*
- *Tip design adds vertical dimension, placing less stress on wound margins*
- *Microbevels provide stability and create a smooth, well-controlled, three-dimensional incision which places less stress on the wound margins, allows the wound to be more flexible during lens insertion with potentially few intra-corneal tears, and helps to promote rapid healing*



Manufactured by



All trademarks are property of their respective owners. ©2016 Surgical Specialties Corporation
All Rights Reserved SSC-098 Rev0 5/19

Distributed by

FCI Ophthalmics

Tel: 800.932.4202

orders@fci-ophthalmics.com

